# ACTIFLUOR™ ADAMTS13 Activity Assay\*

REF 812

EC REP american diagnostica GmbH  
Kaplaneigasse 35, D-64319 Pfungstadt, Germany

## ENGLISH

### INTENDED USE

The ACTIFLUOR™ ADAMTS13 Activity Assay is a fluorescence resonance energy transfer (FRET) assay for the measurement of ADAMTS13 in human plasma. The assay is for *in vitro* diagnostic use.

### EXPLANATION OF THE TEST

ADAMTS13, also known as von Willebrand Factor (vWF) cleaving protease, is a zinc metalloproteinase that cleaves ultra large vWF multimers (UL-vWF) at the Tyr(1605)-Met(1606) bond located in the A2 region of vWF.<sup>1</sup> Studies have shown that low levels of ADAMTS13 activity are associated with Thrombotic Thrombocytopenic Purpura (TTP), a life-threatening hematological condition characterized by low platelet count, microvascular thrombi, red cell fragmentation, CNS and renal complications.<sup>2,3</sup> An ADAMTS13 activity level below 5% of normal leads to an accumulation of UL-vWF multimers in plasma.<sup>4</sup> The UL-vWF multimers bind to receptors on platelets inducing platelet aggregation and formation of intravascular thrombi.

Congenital TTP is a rare inheritable disease caused by mutations within the ADAMTS13 gene which result in the production of non-functional ADAMTS13 protein.<sup>5,6</sup> The acquired form of TTP is an autoimmune-like disorder caused by the development of autoantibodies to ADAMTS13 that inhibit enzyme activity.<sup>7</sup> Quantitation of the ADAMTS13 activity level has been useful in identifying patients with TTP from other thrombocytopenic conditions such as hemolytic uremic syndrome (HUS) and idiopathic thrombocytopenic purpura (ITP).<sup>8</sup>

Assays for measuring ADAMTS13 activity including multimer analysis, collagen binding and vWF ELISAs are long, difficult to perform and highly variable. A rapid and easy to perform FRET assay using a synthetic VWF73 peptidyl substrate was described by Kokame, *et al.*<sup>9</sup> Recently, Zhang, *et al.*, developed a novel recombinant VWF86 FRET substrate based upon autoquenching of fluorescein.<sup>10</sup>

### PRINCIPLE OF THE METHOD

ACTIFLUOR ADAMTS13 is a FRET assay that measures the amount of ADAMTS13 activity in human plasma. A citrated plasma sample is assayed for ADAMTS13 protease activity using recombinant ALEXA488-VWF86 FRET substrate. Proteolytic cleavage of the ALEXA488-VWF86 FRET substrate between the tyr/met residues by ADAMTS13 uncouples the ALEXA fluorochromes resulting in an increase in fluorescence. The increase in fluorescence over time ( $V_{max}$ ) due to cleavage of the substrate by ADAMTS13 is monitored at 37°C using a spectrofluorometer (Ex=485 nm; Em=535 nm). A standard curve is constructed using a normal plasma with a known concentration of ADAMTS13. The ADAMTS13 activity in the plasma is determined by interpolation of the  $V_{max}$  values from the standard curve.

**REAGENTS** – Sufficient for assaying 16 samples in duplicate

- R1 Fluorescence Microwells**, 48 (white)
- R2 ADAMTS13 Standard**, 250 µL, 2 vials (lyophilized)
- R3 ADAMTS13 Positive Control**, 150 µL (lyophilized)
- R4 DMSO**, 0.5 mL
- R5 ALEXA488-VWF86 FRET Substrate**, 300 µL (lyophilized)
- R6 Assay Buffer**, 10 mL, 2 vials (lyophilized)
- R7 ADAMTS13 Inactivated Plasma**, 0.5 mL, 3 vials (lyophilized)

#### WARNING

Source material for some of the reagents in this kit is of human origin. The material has been found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus Type 1 and Type 2 (HIV-1, HIV-2) using FDA approved methods. As no known test method can provide complete assurance that products derived from human blood will not transmit HBsAg, HCV, HIV-1, HIV-2 or other blood-borne pathogens, reagents should be handled as recommended for any potentially infectious human specimen.

#### REAGENT RECONSTITUTION AND STORAGE

##### R2 ADAMTS13 Standard

Add 250 µL of filtered deionized or distilled water to a vial. Mix gently and let stand for 10 minutes at room temperature (18°-25°C). See the vial label for the lot specific standard concentration “C”. Discard unused reconstituted ADAMTS13 Standard after performing the assay.

##### R3 Positive Control

Add 150 µL of filtered deionized or distilled water to the vial. Mix gently and let stand for 10 minutes at room temperature (18°-25°C). The Positive Control is stable for up to one year when stored frozen at –80°C.

##### R5 ALEXA488-VWF86 FRET Substrate

Dilute R4, the DMSO 1:4 with filtered deionized water to prepare a 25% DMSO solution. Add 300 µL of the 25% DMSO solution to the vial. Do not vortex and let stand for 5 minutes. Mix briefly. NOTE: Protect from light. ALEXA488-VWF86 FRET Substrate is stable for up to one year when stored in the dark at –20°C.

##### R6 Assay Buffer

Add 10 mL of distilled water to a vial and let stand for 5 minutes; vortex briefly. A vial contains sufficient Assay Buffer to perform 24 assays. Assay Buffer is stable for up to one month when stored at -20°C.

##### R7 ADAMTS13 Inactivated Plasma

Add 500 µL of distilled water to a vial. Gently mix and let stand for 10 minutes at room temperature (18°-25°C). ADAMTS13 Inactivated Plasma is stable for 4 hours when stored at room temperature or for up to 6 months when stored at –80°C.

#### SPECIMEN COLLECTION AND PREPARATION

Only citrate collected platelet poor plasma may be used for this assay. Do Not Use EDTA. See "Collection, Transport and Processing of Blood Specimens for Testing Plasma-based Coagulation Assays; Approved Guidelines-Fourth Edition", NCCLS Document H21-A4, Vol. 23, No. 35, December 2003. Plasma collection should be performed as follows:

1. Collect 9 parts of blood into 1 part of 3.2% (0.109 M) trisodium citrate anticoagulant solution.
2. Centrifuge the blood sample at 10,000 x g for 15 minutes.
3. Plasma should be stored at room temperature and assayed within 2 hours. Alternatively, plasma may be stored at –80°C for up to 6 months
4. Frozen plasma should be thawed rapidly at 37°C. Thawed plasmas should be stored at room temperature and assayed within 2 hours.

#### PROCEDURE

**Materials Provided** – See Reagents

##### Materials Required But Not Provided

- 0.22 µm filtered deionized H<sub>2</sub>O or distilled water
- 8-channel pipette covering 50-200 µL
- 1-channel pipette covering 10-200 µL
- microfuge tubes, 1.5 mL
- microwell plate fluorometer

#### Preparation of Standards and Reagents

##### A. Preparation of ADAMTS13 Standards

Prepare seven (7) ADAMTS13 standards for generating the standard curve by making serial dilutions of the ADAMTS13 Standard in microfuge tubes as follows:

Tube	Standard Concentration	Volume of	Added to
1	“C” ng/mL	120 µL of Standard	0 µL
2	C/2 ng/mL	60 µL from Tube 1	60 µL of ADAMTS13 Inactivated Plasma
3	C/4 ng/mL	60 µL from Tube 2	60 µL of ADAMTS13 Inactivated Plasma
4	C/8 ng/mL	60 µL from Tube 3	60 µL of ADAMTS13 Inactivated Plasma
5	C/16 ng/mL	60 µL from Tube 4	60 µL of ADAMTS13 Inactivated Plasma
6	C/32 ng/mL	60 µL from Tube 5	60 µL of ADAMTS13 Inactivated Plasma
7	0 ng/mL	0 µL	60 µL of ADAMTS13 Inactivated Plasma

##### B. Preparation of Sample Plasma and Positive Control

Prior to assay, sample plasma and Positive Control are diluted 1: 2 with ADAMTS13 Inactivated Plasma. In a plastic tube, add 50 µL of citrated plasma or reconstituted

Positive Control to 50  $\mu\text{L}$  of ADAMTS13 Inactivated Plasma and mix gently. Keep at room temperature.

### C. Preparation of ALEXA488-VWF86 FRET Substrate


Prewarm Assay Buffer to 37°C prior to performing the assay. Dilute reconstituted Substrate 1:25 with the prewarmed Buffer just before use in Step 5. Note: The assay requires 1.6 mL of prepared substrate for every 16 microwells used.

#### Assay Procedure (kinetic method)

1. Add 80  $\mu\text{L}$  of Assay Buffer to each microwell to be used.
2. Add 20  $\mu\text{L}$  of each ADAMTS13 Standard from Step A to wells A1-A2 through G1-G2 respectively. Add 20  $\mu\text{L}$  of 1:2 diluted Positive Control from Step B to wells H1-H2.
3. Add 20  $\mu\text{L}$  of 1:2 diluted plasma samples from Step B to each of two microwells.
4. Place the microwell plate in a fluorescence plate reader set at Ex=485 nm and Em=535 nm, with a cut-off at 530 nm, and the temperature set at 37°C. For readers using filters, use your FITC filter combination set. Allow the standards and samples to warm for 3 minutes.
5. Add reconstituted ALEXA488-VWF86 Substrate at 1:25 dilution to the prewarmed Assay Buffer from Step C (e.g. add 0.2 mL of Substrate to 4.8 mL of Assay Buffer).
6. Add 100  $\mu\text{L}$  of the 1:25 diluted ALEXA488-VWF86 FRET Substrate to each microwell. Mix plate briefly before making fluorescence readings.
7. Measure the increase in fluorescence for 20 minutes, collecting data at 20 second intervals. Calculate the rate of change in fluorescence, RFU/sec or  $V_{\text{max}}$ .

### RESULTS

A standard curve is constructed by plotting the RFU/sec or  $V_{\text{max}}$  for each standard versus the corresponding concentration of ADAMTS13 in ng/mL. A standard curve should be generated each time the assay is performed. A representative standard curve using a Spectramax® Spectrofluorometer (Molecular Devices, CA) is shown below.

### CALCULATION (Kinetic method)

Interpolate the ADAMTS13 concentration in the assayed sample or Positive Control directly from the standard curve via the RFU/sec or  $V_{\text{max}}$ . Multiply this concentration by 2 to obtain the ADAMTS13 concentration in the sample or Positive Control. If the RFU/sec or  $V_{\text{max}}$  of the sample falls above the range of the standard curve, dilute the plasma further with ADAMTS13 Inactivated Plasma and repeat the assay.


#### Alternative Assay Procedure (Endpoint Method)

If the fluorescence plate reader is not suitable for measurement in the kinetic mode, the endpoint method can be used.

1. Add 80  $\mu\text{L}$  of Assay Buffer to each microwell to be used.
2. Add 20  $\mu\text{L}$  of each ADAMTS13 Standard from Step A to wells A1-A2 through G1-G2 respectively. Add 20  $\mu\text{L}$  of 1:2 diluted Positive Control from Step B to wells H1-H2.
3. Add 20  $\mu\text{L}$  of 1:2 diluted plasma samples from Step B to each of two microwells.
4. Place the microwell plate in a fluorescence plate reader set at Ex=485 nm and Em=535 nm, with a cut-off at 530 nm, and the temperature set at 37°C. For readers using filters, use your FITC filter combination set. Allow the standards and samples to warm for 3 minutes.
5. Add reconstituted ALEXA488-VWF86 Substrate at 1:25 dilution to the prewarmed Assay buffer from Step C (e.g. add 0.2 mL of reconstituted Substrate to 4.8 mL of Assay Buffer).
6. Add 100  $\mu\text{L}$  of the 1:25 diluted ALEXA488-VWF86 FRET Substrate to each microwell. Mix plate briefly before making fluorescence readings.
7. Immediately read the fluorescence at the start of the reaction ( $T=0$ ).
8. Let the plate incubate for 20 minutes at 37°C, mix and read again ( $T=20\text{min}$ ). Calculate the delta RFU ( $\text{delta RFU} = \text{RFU}_{T=20\text{min}} - \text{RFU}_{T=0}$ ).

### RESULTS

A standard curve is constructed by plotting the delta RFU for each standard versus the corresponding concentration of ADAMTS13 in ng/mL. A standard curve should be generated each time the assay is performed. A representative standard curve using a Spectramax® Spectrofluorometer (Molecular Devices, CA) is shown below.

#### CALCULATION (Endpoint method)

Interpolate the ADAMTS13 concentration in the assayed sample directly from the standard curve via the delta RFU of the sample or Positive Control. Multiply this concentration by 2 to obtain the ADAMTS13 concentration in the sample or Positive Control. If the delta RFU of the sample falls above the range of the standard curve, dilute the plasma further with ADAMTS13 Inactivated Plasma and repeat the assay.

#### EXPECTED VALUES

In a study of normal individuals (n=234), the mean ADAMTS13 activity level was found to be  $666 \pm 135$  ng/mL.

#### QUALITY CONTROL

The Positive Control should be included whenever the assay is performed. If the Positive Control fails to be measured within the specified range, the run should be repeated. Contact American Diagnostica if the control repeatedly fails to be measured within the specified range.

#### TRACEABILITY OF CALIBRATORS AND CONTROL MATERIAL

Information on traceability of calibrators and control material is available upon request.

#### PERFORMANCE CHARACTERISTICS

##### Precision

The intra-assay coefficient of variation was determined to be 4.1%

The inter-assay coefficient of variation was determined to be 4.4%.

#### REFERENCES

1. Furlan, M., Robles, R. and Lämmle, B. L. *Blood* 1996, **87**: 4223-4234.
2. Tsai, H. M. *Blood* 1996, **87**: 4235-4244.
3. Furlan, M., *et al.* *N Engl J Med* 1998, **339**: 1578-1584.
4. Moake, J. L., *et al.* *N. Engl J Med* 1982, **307**: 1432-1435.
5. Zheng X, *et al.* *J Biol Chem* 2001, **276**: 41059-41063.
6. Levy, G. G., Nichols, W. C. and Lian, E. C. *Nature* 2001, **413**: 488-494.
7. Tsai, H. M., Rice, L. and Sarode, R. *Ann Intern Med* 2000, **132**: 794-799.
8. Unpublished data from American Diagnostica Inc.
9. Kokame K, Nobe Y, Kokubo Y, Okayama, A. and Miyata, T. *Br. J. Haematol.* 2005, **129**: 93.
10. Zhang, L., Lawson, H. L., Harrish, V. C., Huff, J. D., Knovich, M. A. and Owen, J. *Anal Biochem* 2006:358:298-300.

Spectramax is a trademark of Molecular Devices, Sunnyvale, CA, USA

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8. Unpublizierte Daten von American Diagnostica Inc.