american diagnostica inc.

500 West Avenue, Stamford, CT 06902
Tel. (203) 602-7777 • Fax. (203) 602-2221

ENGLISH

ACTICHROME[®] Heparin (anti-fIIa)

REF 820

(CPT Code No. 85520)

american diagnostica GmbH
Kaplaneigasse 35, D-64319 Pfungstadt, Germany

INTENDED USE

ACTICHROME[®] Heparin (anti-fIIa) is an amidolytic chromogenic assay intended for the quantitative determination of therapeutic heparin in human plasma by measurement of factor IIa (thrombin) activity.

PRINCIPLE

The inhibitory effect of antithrombin III (AT-III) on factor IIa (thrombin), factor Xa, and other coagulation serine proteases in plasma is increased several thousand-fold by heparin. This inhibition accounts for the anticoagulant effect of heparin. The quantitative determination of plasma heparin levels by the measurement of their anti-factor IIa activity is a necessary tool for monitoring treatment efficacy.

Low molecular weight heparin (LMWH) preparations appear to catalyze the reaction between thrombin (fIIa) and AT-III less readily than the reaction between factor Xa and AT-III. Unfractionated heparin (UFH) catalyzes both reactions equally. In this assay, the rate of factor IIa inhibition is directly proportional to the heparin concentration since both factor IIa and AT-III are in excess. The residual factor IIa activity is inversely proportional to the heparin concentration.^{1,2}

REAGENTS

The kit contains sufficient reagents to perform 48 tests using a manual end-point method. Automated instruments yield a higher number of tests.

- R1 **Bovine Thrombin Reagent:** 6 vials each containing bovine thrombin (lyophilized).
- R2 **Human Antithrombin III Reagent:** 6 vials each containing human antithrombin III (lyophilized).
- R3 **SPECTROZYME[®] TH:** 6 vials each containing 1.8 µmoles of substrate (lyophilized).

REAGENTS AND EQUIPMENT REQUIRED BUT NOT PROVIDED

spectrophotometer set at 405 nm and laboratory timer or,
automated chromogenic photo-optical instrument
0.22 µm filtered deionized/distilled water for dilution
0.9% NaCl
glacial acetic acid
25 µL, 200 µL, 500 µL and 2500 µL pipettes
37°C water or dry bath
plastic test tubes or cuvettes
linear graph paper
commercially available normal plasma (e.g. ADI Cat. No. 258N)

REAGENT RECONSTITUTION AND STABILITY

Unreconstituted reagents are stable until the expiration date indicated on the label when stored at +2° - +8°C.

R1 Bovine Thrombin Reagent: Reconstitute with 2 mL filtered deionized/distilled water. Reconstituted reagent is stable for 2 weeks at +2° - +8°C and for 4 months at -20°C.

R2 Human Antithrombin III Reagent: Reconstitute with 5 mL filtered deionized/distilled water. Reconstituted material is stable for 2 weeks at +2° - +8°C and for 4 months at -20°C.

R3 SPECTROZYME® TH: Reconstitute with 2 mL filtered deionized/distilled water. Reconstituted material is stable for 2 months at +2° - +8°C and for 6 months at -20°C.

WARNING

The source material for the human antithrombin III has been found to be non-reactive for Hepatitis B Surface Antigen (HBsAg), Hepatitis C Virus (HCV) and Human Immunodeficiency Virus Type 1 and Type 2 (HIV-1, HIV-2) using FDA Approved methods. As no known test method can provide complete assurance that products derived from human blood will not transmit HBsAg, HCV, HIV-1, HIV-2 or other blood-borne pathogens, this reagent should be handled as recommended for any potentially infectious human specimen.

SPECIMEN COLLECTION

Nine volumes of blood is collected in 1 volume of 0.1M trisodium citrate and centrifuged at 3,000 x g for 15 minutes (See NCCLS standard H21-A3). Exercise great care in the centrifugation of the samples in order to minimize platelet alteration, as well as to achieve maximum removal of the platelets. Perform centrifugation within one hour after sample collection. Plasma should be stored at +2° - +8°C and assayed within 2 hours. Platelet poor plasma may be stored at -20°C for up to one month and then thawed once at 37°C, 30 minutes before testing.

ASSAY CALIBRATION

Pooled normal human plasma (from at least 20 normal donors) that has been collected in the same manner as plasmas to be tested may be used for preparation of the heparin standards. The specific heparin preparation used during therapy must be added to the pooled normal plasma to prepare the standards.

Prepare a 6.0 USP unit/mL solution of heparin in saline (0.9% NaCl). When starting with a 1000 unit/mL heparin stock solution, pipette 30 µL of the heparin stock into 4.97 mL of the saline. Then prepare the heparin standards using the 6.0 USP unit/mL solution as shown in Table 1.

Table 1

Heparin Standard	Pooled Normal Plasma	Volume of Heparin Standard
0.6 unit/mL	900 µL	100 µL of 6.0 USP unit/mL
0.3 unit/mL	500 µL	500 µL of 0.6 USP unit/mL
0.15 unit/mL	500 µL	500 µL of 0.3 USP unit/mL
0.0 unit/mL	500 µL	None

Prepare a 1:16 dilution of each patient plasma specimen and each heparin standard by adding 25 µL of plasma specimen or standard to 375 µL of Human Antithrombin III Reagent.

PROCEDURE

Manual End-point Method

1. Add 200 µL of Human Antithrombin III Reagent to a plastic test tube.
2. Add 25 µL of 1:16 dilution of plasma sample or heparin standard.
3. Mix and incubate at 37°C for 2 minutes.
4. Add 200 µL of Bovine Thrombin Reagent.
5. Mix and incubate at 37°C for exactly 2 minute.
6. Add 200 µL of SPECTROZYME TH Substrate.
7. Mix and incubate at 37°C for exactly 1 minute.
8. Add 200 µL of glacial acetic acid. Mix.
9. Add 200 µL of filtered deionized water (optional)*.
10. Read absorbance at 405 nm in a 1 cm semi-micro cuvette against a blank prepared in the following order:
 - 200 µL acetic acid
 - 200 µL Human Antithrombin III Reagent
 - 25 µL pooled normal plasma diluted 1:16 in Human Antithrombin III Reagent
 - 200 µL Bovine Thrombin Reagent
 - 200 µL SPECTROZYME TH Substrate
 - 200 µL water (optional)*


* Some spectrophotometers require a minimum of 1 mL volume in cuvette.

Automated Instrumentation Method

Automated chromogenic photo-optical instrument procedures utilize an initial rate reaction rather than an end-point method as described above. Please refer to the manufacturer's instrument manuals for chromogenic assay procedures. The performance of this assay was determined on an MLA® 900C optical analyzer. The instrument application is available from American Diagnostica upon request.

CALIBRATION CURVE

Plot the absorbance obtained at 405 nm for each heparin standard against its concentration on linear graph paper. The heparin in patient plasma can be determined by interpolation from the calibration curve (see sample calibration curve). The sample curve shown below is for example only. A standard curve should be generated each time the assay is performed.

LIMITATIONS OF THE PROCEDURE

This assay is not optimized to measure patient samples containing low molecular weight heparins (LMWH). LMWH preparations appear to catalyze the reaction between thrombin (fIIa) and AT-11 less readily than the reaction between factor Xa and AT-III.

Different heparin preparations may elicit different slopes in the standard curve. The heparin administered to the patient must be used to generate the standard curve.

Platelet contamination in plasma samples may result in the release of platelet factor 4 (PF4), a potent heparin neutralizer. PF4 released into plasma may result in an underestimation of the heparin concentration. Freezing and thawing of platelets is known to release PF4 from the platelets. If plasma is to be frozen and thawed prior to use in the assay, then the plasma must be prepared as described in the sample preparation section so that it is platelet poor.

INTERFERENCES

Icteric, lipemic and hemolyzed samples may interfere with assay results. Studies have shown that these samples result in an underestimation of the heparin concentration.³

PERFORMANCE CHARACTERISTICS

Accuracy

A study performed using an automated instrument to evaluate ACTICHROME Heparin (anti-fIIa) against another commercially available chromogenic heparin assay gave a correlation coefficient of 0.967 for 88 samples tested, with a regression equation of $y = 0.845x + 0.027$. The Standard Error of Estimation = 0.03.

Precision

The following estimates of coefficient of variation (CV) were observed using an automated chromogenic photo-optical instrument.

Table 2

Heparin Concentration	Intra-Assay Variation	Inter-Assay Variation
0.42 units/mL	4.7%	7.6%
0.23 units/mL	10.8%	9.6%

Sensitivity

ACTICHROME Heparin (anti-fIIa) is designed to give a linear standard curve for heparin levels between 0.0 and 0.6 units/mL.

Specificity

Specificity of the assay has been ensured by the use of purified bovine thrombin and purified human AT-III. In addition, the chromogenic substrate used to measure residual thrombin activity is highly specific for thrombin.

Quality Control

Quality control samples should be run according to Laboratory Guidelines. Control samples should be run on freshly reconstituted reagent vials, and each time the assay is run. Commercial heparin control plasma may be used for quality control of the assay. Alternatively, pooled normal plasma containing 0.5 and 0.25 units/mL heparin can be made and used as in-house controls. If quality control samples fail to yield heparin concentrations within the correct ranges in any given test, the test should be repeated with freshly reconstituted reagents.

REFERENCES

1. Teien, A. N., Lie, M. and Abildgaard, U. Assay of heparin in plasma using a chromogenic substrate. *Thrombosis Research* 1976, 8: 413-416.
2. Teien, A. N. and Lie, M. Evaluation of an amidolytic heparin assay method: Increased sensitivity by adding purified antithrombin III. *Thrombosis Research* 1977, 10: 399-410.
3. Data on file at ADI.

MLA is a registered trademark of Instrumentation Laboratory, SpA

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