



american diagnostica inc.

## sheep anti-human factor VII IgG

Product No. 231S

Lot No. 080314

Exp. 2013-03-14

**Preparation:** Ammonium sulfate precipitation followed by ion-exchange chromatography and gel-filtration

**Specific Activity:** Not determined.

**Extinction Coefficient:**  $\Sigma_{1\%, 1 \text{ cm}, 280 \text{ nm}} = 14.0$

**Molecular Weight:** 150,000 Daltons

**Concentration:** 4.6 mg/mL

**Presentation:** 1 mg in 217.4  $\mu\text{L}$  of a 50% glycerol/ $\text{H}_2\text{O}$  solution, vialled in a plastic screw-capped cryotube.

**Storage:** Store at  $-20^\circ\text{C}$