



american diagnostica inc.  
*Creating products for life<sup>®</sup>*

## sheep anti-human factor XIII IgG

### Product No. 4130

#### Immunogen/ Description

Antisera from sheep immunized with purified intact human plasma factor XIII.

#### Preparation

The IgG fraction purified from crude anti-sera via ammonium sulphate precipitation and gel filtration.

#### Properties

Reacts with the a and b subunits of human factor XIII per immunoblot studies.

Is useful in both ELISA and Western blot protocols.

May be used for immunopurification of human factor XIII.

#### Presentation

Screw-capped glass vial containing 1 mg of purified IgG lyophilized from 1 mL of a 0.15M PBS, 1% mannitol, pH 7.4 buffer.

#### Reconstitution

Add 1 mL of filtered deionized water for a 1 mg/mL stock solution.

#### Storage

Store lyophilized IgG at 2° - 8°C. Reconstituted IgG should be aliquoted and stored at -20°C or colder.

#### Related Products

Human Factor XIII, Cat. No. 413. Human Factor XIIIa, Cat. No. 413a.