



american diagnostica inc.  
*Creating products for life<sup>®</sup>*

## **goat anti-human factor IX IgG**

**Product No. 499   Lot No. 030327**

### **Immunogen**

Goat immunized using purified human plasma factor IX. Antisera pool prepared from different goats and bleeds.

### **Purification**

Crude antiserum purified using immobilized Protein A, to isolate the IgG fraction.

### **Properties**

Monospecific antiserum recognizing zymogen and protease forms of human plasma factor IX, in addition to a commercially available recombinant preparation, per Ouchterlony studies. Animal plasma fIX less strongly recognized

### **Applications**

Ouchterlony studies of factor IX preparations, counter and bidimensional electrophoresis, nephelometry and immunopurification.

### **Presentation**

Plastic cryo-tubes containing 1 mg of purified goat IgG in 78  $\mu$ L of a solution of 50% glycerol/10 mM Hepes, 150 mM NaCl, pH 7.4. 15M PBS, pH 7.4

### **Storage**

Store the vial at -20°C. (Ambient shipping conditions are acceptable.)

### **Related Products**

Various different monoclonal antibodies against human factor IX. Bovine factor IX, Human factor IX zymogen and protease.

© ADI990113.rev2.030327