



**HRP labeled monoclonal antibody  
against human factor VIII:C**

**Product No. ESH-8PO**

**American Diagnostica GmbH**

Kaplaneigasse 35 • D-64319 Pfungstadt • Germany  
Tel: +49 (0)6157 990899 • Fax: +49 (0)6157 990808

## Description

HRP (horseradish peroxidase) labeled monoclonal antibody ESH-8 directed against human Factor VIII:C. ESH-8 is reactive with the C2 domain of the FVIII light chain (aa: 2248-2285).

## Applications

Useful as a secondary antibody in an ELISA format.

## Presentation

Screw-capped plastic vial containing 100 µg of HRP-labeled antibody in 200 µL stabilizing buffer. Concentration approx. 0.5 mg/ml.

## Storage

Store the antibody at 2 – 8 °C, do not freeze.

## References

1. Griffin, B. D., *et al.* *Thrombosis and Haemostasis*, 1986; **55**: 40-46.
2. Hornsey, V. S., *et al.* *Transfusion Medicine*, 1992; **2**: 223-229.
3. Zatloukal, K., *et al.* *Proceeds of the National Academy of Science USA*, 1994; **91**: 5148-5152.
4. Scandella, D., *et al.* *Blood*, 1995; **86**: 1811-1819.

## Related Products

Native and Recombinant human factor VIII. Native vWF, vWF mono and polyclonal antisera and the IMUBIND® vWF ELISA (No. 828), IMUBIND® vWF Activity ELISA (No. 885), and the IMUBIND® fVIII ELISA (No. 884CON).

©ADG12052011